

Extend your range.

Motion Foot® SLX



**Foot
Weight**
23.8 oz.
675 g

**Maximum
Patient Weight**
276 lbs.
125 kg

**Build
Height**
4.1 in.
10.4 cm

**Activity
Level**
K3

Fillauer®





Hydraulic Motion with Less Weight

The Motion Foot SLX allows for 13° of plantarflexion and 3° of dorsiflexion in a foot that's lightweight, low build height, and more durable than ever, providing compliance and confidence on slopes and uneven terrain.

Stable & Smooth

The range of motion provides stability in early stance, particularly for transfemoral users, and the dorsiflexion stop allows for energy storage and initiation of a ground reaction force that stabilizes the knee from mid-stance to terminal stance.



Dialed-in Response

The SLX hydraulic unit has individual plantarflexion and dorsiflexion adjustments for smooth motion throughout the gait cycle.

Motion Foot® SLX

Specifications

Patient Weight	275 lbs. (125 kg)
Foot Build Height	4.08 in. (10.36 cm)
Foot Weight (26 cm, Category C)	23.8 oz. (675 g)
Functional Level	K3
Primary Materials	Carbon composite, stainless steel, titanium, and aluminum
Durable	Meets ISO-10328 Standard
Waterproof	The foot unit is waterproof to 1 meter.
Warranty	36 months



Ordering Information

Foot includes Spectra Sock. Foot Shell sold separately.

213-10-XX-XX Motion Foot SLX foot, foot sizes 22 – 30 cm

Patient Weight		22 cm	23 cm	24 cm	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm
100 – 130 lbs. 45 – 59 kg		22-A1	23-A1	24-A1	25-A3	26-A3				
131 – 161 lbs. 60 – 73 kg		22-A1	23-A1	24-A1	25-B3	26-B3	27-C5	28-C5		
162 – 192 lbs. 74 – 87 kg		22-B2	23-B2	24-B2	25-C4	26-C4	27-D5	28-D5	29-E6	30-E6
193 – 225 lbs. 88 – 102 kg		22-B2	23-B2	24-B2	25-D4	26-D4	27-E6	28-E6	29-F7	30-F7
226 – 260 lbs. 103 – 118 kg					25-E4	26-E4	27-F6	28-F6	29-G7	30-G7
261 – 276 lbs. 119 – 125 kg							27-G7	28-G7	29-G7	30-G7

Foot Shell

Made with BIOVYN™, a bio-attributed PVC that is based on 100% renewable feedstock, our product not only contributes to significantly reduced carbon emissions but also offers enhanced durability. To complement the foot shell, a sleek, black prosthetic sock, sized appropriately, is now included.



Specifications

Product Weight (26 cm)	Material	Warranty
8.9 oz. (251 g)	BIOVYN™ PVC	9 months

Ordering Information

Replace **XX** with foot length (22 – 30 cm). Foot shell includes black prosthetic sock.

Color	Left	Right
03	45-XX-13-03-3	46-XX-13-03-3
09	45-XX-13-09-3	46-XX-13-09-3
13	45-XX-13-13-3	46-XX-13-13-3
99	45-XX-13-99-3	46-XX-13-99-3

Prosthetic Socks

Size	Part Number
22 – 24 cm	636202224
25 – 27 cm	636202527
28 – 30 cm	636202830

Fillauer®

Fillauer LLC

2710 Amnicola Highway
Chattanooga, TN 37406
423.624.0946