



# TASKA<sup>®</sup>





The world's only rugged,  
water-resistant, multi-  
articulating hand.



## FLEXIBLE FINGERS

The fingers on the Taska hand are flexible and can laterally spread to allow the hand to effectively grip a wide variety of objects firmly.

## CONFIDENT GRIP

Hold things with more confidence. Soft pads between the fingers and on fingertips, complement the flexible fingers to provide unparalleled grip security.

## MOTORIZED THUMB ROTATION

High speed thumb rotation mimics the motion of the natural hand and allows nearly instant change of grip pattern.

## THE TASKA® LOW PROFILE WRIST OPTION

A low profile wrist option extends water resistance up the forearm. The LP Wrist also allows for 90° of passive rotation.

## KNUCKLE BREAKAWAY


Knuckles break away under large loads to prevent damage and are user resettable to provide additional robustness.

## GRIP CYCLE BUTTONS

Simple touch-and-go, easy-to-use grip cycle buttons provide access to grips. EMG grip patterns can also be disabled at the touch of a button.

## INTEGRATED FLEXIBLE WRIST

Every model of the Taska hand has a built-in, flexible wrist with three easily lockable positions, as well as center-biased free flexion.



8 ¼



Black



Sand



White

Motion Control  
division of  
*Fillauer*®



7 ¾

801.326.3434

Exclusive Distributor for U.S and Canada

[www.fillauer.com/taska](http://www.fillauer.com/taska)